

Multi-centre Study of Doxorubicin in the **Treatment of Hepatocellular Carcinoma** by DEBDOX™ with DC Bead™

PRECISION V Trial Design

- Control arm of conventional transarterial chemoembolisation (cTACE) with doxorubicin
- 200 patients (100 per arm)
- 23 European centres in Austria, France, Germany, Greece, and Switzerland
- Patients receive up to 3 treatments (baseline, 2 and 4 months)
- Follow-up period: 6 months

Patient Demographics

Characteristics	DC Bead® (n=102)	cTACE (n=110)	
Age Mean (±sd)	67.3 years (±9.1)	67.4 years (±8.8)	
Sex (Male/Female)	79/14	95/13	
Aetiology (HepC/HepB/Alcohol/Other)	22/16/43/21	18/18/57/25	
Okuda (I/II)	79/14	103/5	
BCLC (A/B/C)*	24/69/0	29/79/0	
No. Lesions Median (range)	2.8 (1-20)	3.8 (1-50)	
Sum Longest Diameter Mean (±sd)	88.9mm (±52.1)	89.2mm (±59.3)	
Liver Involvement Mean (range)	16.1% (<10-50)	16.1% (<10-50)	

*BCLC Classification according to tumour stage (Llovet et al Lancet 2003)

Product, Dose and Technique Guidelines

DEBDOX™ with DC Bead™

- 2 x 2ml vials of DC Bead™ (total 4ml) loaded at 37.5mg/ml for total dose of 150mg
- 1 vial of 300-500μm followed by 1 vial of 500-700μm

cTACE

- Doxorubicin dose of 50-75mg/m² to maximum of 150mg
- Physician preference for embolic

Technique for both groups

- Unifocal tumours treated with selective segmental chemoembolisation
- Microcatheter could be used
- Bilobar disease: both lobes treated within a 3-week period
- Embolisation to stasis in 2nd or 3rd order branches
- DC Bead[™] group: additional Bead Block[™] could be used

DCBead



DCBead

DEBDOX[™] for Hepatocellular Carcinoma

DC Bead™ DEBDOX™ Bibliography

emoembolization With Doxorubicin-Eluting Beads for Unresectable Hepatocellular Carcinoma: e-Year Survival Analysis gari K, Pomoni M, Moschouris H et al. Cardiovasc Intervent Radiol (2012) DOI: 10.1007/s00270-012-0394-0 urvival of Patients with Hepatocellular Carcinoma Treated by Transarterial Chemoembolization (ACE) using DC Beads. Implications for Clinical Practice and Trial Design rel M, Reig M, Forner A et al. J Hepatol 56 (2012): 1330-1335 ranscatheter Arterial Chemoembolization for Liver Cancer: Is It Time to istinguish Conventional from Drug-Eluting Chemoembolization? api E and Geschwind JF, Cardiovasc Intervent Radiol 31 (2010); Published online nparison of Conventional Transarterial Chemoembolization (TACE) and Chemoembolization Doxorubicin Drug Eluting Beads (DEB) for Unresectable Hepatocellular Carcinoma sekaran R, Kooby D, Staley C et al. J of Surg Onc 101 (2010): 476-480 ontrast-Enhanced Ultrasonography of Hepatocellular Carcinoma After Chemoembolisation sing Drug-Eluting Beads: A Pilot Study Focused on Sustained Tumor Necrosis houris H. Malagari K. Papadaki MG et al. Cardiovasc Intervent Radiol 33 (2010): 1022-1027 Prognostic Factors for Survival in Patients with Unresectable Hepatocellular Carcinoma Undergoing Chemoembolization with Doxorubicin Drug-Eluting Beads: a Preliminary Study aran R, Kooby D, Staley C et al. HPB Oxford 12 (2010): 174-180 co-regional Interventional Treatment of Hepatocellular Carcinoma: hniques, Outcomes, and Future Prospects ncioni R. J Hepatol 52 (2010): 762-773 Transarterial Chemoembolization with Epirubicin-eluting Beads versus Transarterial Embolization before Liver Transplantation for Hepatocellular Carcinoma colini A, Martinetti L, Crespi S et al. J Vasc Interv Radiol 21 (2010): 327-332 spective Randomized Comparison of Chemoembolization with Doxorubicin-Eluting ds and Bland Embolization with Bead Block" for Hepatocellular Carcinoma agari K, Pomoni M, Kelekis A et al. Cardiovasc Intervent Radiol 33 (2010): 541-551

DC Bead™ Important Information

Indication: DC Bead™

• DC Bead is primarily intended as an embolic agent for the treatment of malignant hypervascularised tumour(s) DC Bead is compatible with doxorubicin, which can be loaded prior to embolisation and then, as a secondary

This product and/or all indications may not be available in your territory DC Bead is not cleared by the FDA for sale or distribution in the USA.

Cautions: DC Bead™

- Embolisation with DC Bead should only be performed by a physician with appropriate interventional occlusion.
- raining in the region intended to be embolised . Do not use if the vial or packaging appear damaged
- Ensure that DC Bead is an appropriate size for the intended vasculature
- . Consider upsizing to a larger size of DC Bead in the presence of AV shunts or if angiographic evidence of embolisation does not appear quickly during delivery
- Consideration should be given to Tc99m-MAA scanning if there is suspicion of AV shunting
- Exceeding a loading dose of 37.5mg doxorubicin per 1ml DC Bead may lead to some systemic distribution of doxorubicin and related side effects

Potential Complications: DC Bead™

- Undesirable reflux or passage of DC Bead into normal arteries adjacent to the
- targeted lesion or through the lesion into other arteries or arterial beds
- Pulmonary embolisation
- Ischaemia at an undesirable location . Capillary bed saturation and tissue damage
- Ischaemic stroke or ischaemic infarction · Vessel or lesion rupture and haemorrhage
- Neurological deficits including cranial nerve palsies

- · Foreign body reactions necessitating medical intervention · Infection necessitating medical intervention
- $\bullet \ \, \text{Clot formation at the tip of the catheter and subsequent dislodgement causing arterial thromboembolic sequelae} \\$

is a trademark of BTG International Ltd. BTG and the BTG roundel logo are registered trademarks of BTG International Ltd. © Copyright 2009-14 Biocompatibles UK Ltd. EC09-020 Rev 6.

spective Randomized Study of Doxorubicin-Eluting-Bead Embolization in the atment of Hepatocellular Carcinoma: Results of the PRECISION V Study mer J, Malagari K, Vogl T et al. Cardiovasc Intervent Radiol 33 (2010): 41-52

Single-Center Phase II Trial of Transarterial Chemoembolization with Drug-Eluting Beads for Patients with Unresectable Hepatocellular Carcinoma. Initial Experience in the United States yes D, Vossen J, Geschwind J et al. The Cancer Journal 15 (2009): 526-532

Drug-Eluting Bead Therapy in Primary and Metastatic Disease of the Liver ter S and Martin RCG. HPB Oxford 11 (2009): 541-550

Transcatheter Chemoembolization in the Treatment of HCC in Patients not Eligible for Curative Treatments: Mid-term Results of Doxorubicin-loaded DC Bead ladari K. Alexopoulou E, Chatzimichail K et al. Abdominal Imaging 33 (2008): 512-9

oxorubicin-eluting Bead-enhanced Radiofrequency Ablation of Hepatocellular arcinoma: A Pilot Clinical Study cioni R, Crocetti L, Petruzzi P et al. J Hepatol 49 (2008): 217-222

ransarterial Chemoembolization of Unresectable Hepatocellular Carcinoma vith Drug Eluting Beads: Results of an Open-Label Study of 62 Patients agari K, Chatzimichael K, Alexopoulou E et al. Cardiovasc intervent Radiol 31 (2008): 269-280

nemoembolization (TACE) of Unresectable Intrahepatic Cholangiocarcinoma th Slow-Release Doxorubicin-Eluting Beads: Preliminary Results liberti C, Benea G, Tilli M et al. Cardiovasc Intervent Radiol 31 (2008): 883-888

A Phase I/II Trial of Chemoembolization for Hepatocellular Carcinoma Using a Novel Intra-Arterial Drug-Eluting Bead on RTP, Tso WK, Pang RWC et al. Gastroenterol Hepatol 5 (2007): 1100-1108

emoembolization of Hepatocellular Carcinoma with Drug-Eluting Beads: icacy and Doxorubicin Pharmacokinetics rela M, Real MI, Burrel M et al. J Hepatol 46 (2007): 474-481

DC Bead™ Ordering Information:

Label Colour and Size	100-300μm	300-500μm	500-700μm	
Volume of Beads	2ml	2ml	2ml	
Product Code	DC2V103	DC2V305	DC2V507	


Biocompatibles UK Ltd, a BTG International group company Lakeview

Riverside Way Watchmoor Park GU15 3YL, UK

Tel: +44(0) 1276 902 020 email: marketing@btgplc.com www.btg-im.com

Greater tolerability Increased efficacy

for more patients





There are no limits to where our ideas will take us. DC Bead™ and Bead Block™ are manufactured by Biocompatibles UK Ltd, a BTG International group company, DC Bead, Bead Block and DEBDOX are trademarks of Biocompatibles UK Ltd. DC Bead and Bead Block are registered trademarks in the EU and certain other territories. 'Imagine where we can go'

Imagine where we can go.



There are no limits to where our ideas will take us.

Imagine where we can go.

Greater tolerability Increased efficacy

for more patients



Treatment of Hepatocellular Carcinoma

by DEBDOX™ with DC Bead™

"We believe that these results show that DC Bead™ is a better treatment than conventional TACE. An improved response with significantly lower toxicity is unusual for a new cancer therapy."

PRECISION V Publication Committee

Professor Johannes Lammer, Vienna Professor Katarina Malagari, Athens Professor Alban Denys, Lausanne

Professor Riccardo Lencioni, Pisa Professor Frank Pilleul, Lvon Professor Anthony Watkinson, Exeter Professor Thomas Vogl, Frankfurt

PRECISION V Conclusions⁵

- DEBDOX™ with DC Bead™ is well tolerated, efficacious and reproducible
- There is a significant reduction in liver toxicity with DEBDOX™ with DC Bead™
- There is a significant advantage of DEBDOX™ with DC Bead™ in more advanced patients those with more compromised liver function, poorer performance status, bilobar disease and recurrent disease - greater response, greater disease control and improved safety
- Currently AASLD guidelines do not recommend chemoembolisation for Child B and ECOG 1 patients. The PRECISION V data show that these patients can now be safely treated with **DEBDOX™** with DC Bead™

DC Bead™ is an embolic Drug-Eluting Bead capable of loading and releasing in a controlled manner high doses of chemotherapeutic agents.


DC Bead[™] Indication for Use

DC Bead™ is CE-marked and is intended to be loaded with doxorubicin for the purpose of:


- Embolisation of vessels supplying malignant hypervascularised tumour(s)
- Delivery of a local, controlled, sustained dose of doxorubicin to the tumour(s)
- Doxorubicin maximum dose of 37.5mg/ml and 150mg per treatment with 4ml DC Bead™

Interaction of Doxorubicin With DC Bead™ Sulphonate Groups

Hvdrated Beads




Loaded Beads



DC Bead[™] Presentation


- Novel N-fil technology sulphonate modified hydrogel polymer
- Blue tinted to aid visualisation
- Delivered as vials containing 2ml Beads in 6ml saline
- Precise calibration to achieve an accurate level of embolisation







Varela et al

exposure to doxorubicin.






In the PRECISION V clinical trial, patients who received DEBDOX™ with DC Bead™ had a significant (p<0.001) reduction in doxorubicin-related systemic adverse events, despite receiving 30% more doxorubicin.



DC Bead™ improved response in all treated patients. More Advanced Patients demonstrated a significant improvement (p<0.05). DC Bead™ patients experienced fewer treatment-related side effects compared to the control group.

	DC Bead™ Response (%)			cTACE Response (%)			
Classification	DC	OR	CR	DC	OR	CR	
Child Pugh B	63	44	25	32	21	16	
ECOG 1	63	63	37	32	29	14	
Bilobar	59	49	17	49	40	13	
Recurrent Disease	73	55	27	54	31	15	
DC = Disease Control OR = Objective Response CR = Complete Response							



The elevation of liver enzyme levels after each

of the three treatments was significantly less in

patients receiving DC Bead™, demonstrating

that DEBDOX™ with DC Bead™ is less toxic to healthy liver.

^{2.} Varela, M., Real, M.I., Burrel, M. et al. Chemoembolization of hepatocellular carcinoma with drug eluting beads: Efficacy and doxorubicin pharmacokinetics. Journal of Hepatology 2007; 46:474-481.
3. Poon, R.T.P., Tso, W.K., Pang. R.W.C. et al. A phase VII Trial of chemoembolization for hepatocellular carcinoma using a novel intra-arterial drug-eluting bead. Clinical Gastroenterology and Hepatology 2007; 5:100-1108 4. Raoul, J.L., Herebach, D., Bretagne, J.F. et al. Chemoembolisation of hepatocellular carcinomas: a study of the biodistribution and pharmacokinetics of doxorubicin. Cancer 1992; 70:585-590

^{5.} Lammer J, Malagari K, Vogl T et al. Prospective Randomized Study of Doxorubicin-Eluting-Bead Embolization in the Treatment of Hepatocellular Carcinoma: Results of the PRECISION V Study.

^{6.} PRECISION V Clinical Study Final Results. Data on file at BTG International Ltd